



CONTENTS PAGE





ADAPTION SERVICE


















Tenby



Mod Shoe



CHOOSE STANDARD OR LOW OPENING AT NO EXTRA COST ON MOD SHOE!!!

LADIES | SANDALS





LADIES BOOTS



LADIES | BOOTS



LADIES | BOOTS



CHOOSE STANDARD OR LOW OPENING AT NO EXTRA COST

ON MOD BOOT!!!











OUR HISTORY

Ken Hall himself hailed from a long line of shoemakers, prominent in Northamptonshire from the turn of the last century. His main interest was bespoke footwear and, having made the transition from manufacturing for the high street to the manufacture of orthopaedic footwear in the 1970's, he rapidly gained his reputation for producing practical and functional boots and shoes for people with problem feet.


Hall's staff of four soon expanded and in the mid seventies some of the first off the shelf or modular orthopaedic footwear were developed. The company expanded steadily during the eighties and were producing ranges that included adults, children's footwear, plus a full made – to – measure range.

In 1991, bigger premises were required and we transferred to our current factory, which was constructed as a purpose built shoe factory in 1873. This move allowed for improved production facilities and greater stock holding, thus improving delivery times and service, which we maintain to this day.










Regent















Dalton



Purbeck






MENS BOOTS



MENS | BOOTS





MENS | BOOTS



CHOOSE STANDARD OR LOW OPENING AT NO EXTRA COST ON MOD BOOT!!!















Ashby

Effortlessly cool minimalist trainer design that can work well with every outfit. Styled with a hidden padded collar and rubber foxing strip all round as standard. Shown in Black, White and Tan 304 Softy. Lined fully with Microfibre.






SHEARLING LINED HOUSE SHOES & BOOTS



UNLINED FELT SHOE & BOOTS





FLEXIBLE NEOPRENE

Neoprene can be incorporated into most of our designs to allow maximum accommodation.

PERFECT FOR THAT EXTRA FLEXIBILITY

Neoprene material has a good water resistance and maintains flexibility when subject to a wide range of temperatures. Therefore, it is able to resist degradation from ozone, sunlight, oxidation, sand, snow, and dust – essentially all weather conditions.

Below are a few sample ideas:





SIZING CHART | LADIES

Carewear Standard & Carewear Pro-Fit Lasts Narrow Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	78mm	177 – 201mm	170 – 194mm	21-33mm	52
1 1/2	79mm	180 - 204mm	173 - 197mm	21-33mm	53
2	80mm	183 – 207mm	176 – 200mm	22-34mm	54
2 1/2	81mm	186 - 210mm	179 - 203mm	22-34mm	54
3	82mm	189 – 213mm	182 – 206mm	23-35mm	55
3 1/2	83mm	193 - 217mm	186 - 210mm	23-35mm	56
4	84mm	196 – 220mm	189 – 213mm	24-36mm	57
4 1/2	85mm	199 - 223mm	192 - 216mm	25-37mm	57
5	86mm	203 – 227mm	196 – 220mm	25-37mm	58
5 1/2	87mm	206 - 230mm	199 - 223mm	26-38mm	59
6	88mm	210 – 234mm	202 – 226mm	26-38mm	60
6 1/2	89mm	213 - 237mm	205 - 229mm	27-39mm	60
7	90mm	216 – 240mm	208 – 232mm	27-39mm	61
7 1/2	91mm	219 - 243mm	211 - 235mm	28-40mm	62
8	92mm	222 – 246mm	214 – 238mm	28-40mm	63
8 1/2	93mm	225 - 249mm	217 - 241mm	29-41mm	63
9	94mm	229 – 253mm	221 – 245mm	29-41mm	64
9 1/2	95mm	232 - 256mm	224 - 248mm	30-42mm	65
10	96mm	235 – 259mm	227 – 251mm	30-42mm	66

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	83mm	194 – 218mm	188 – 212mm	24-36mm	53
1 1/2	84mm	198 - 222mm	191 - 215mm	24-36mm	54
2	85mm	201 – 225mm	194 – 218mm	25-37mm	55
2 1/2	86mm	204 - 228mm	198 - 222mm	25-37mm	55
3	87mm	208 – 232mm	201 – 225mm	26-38mm	56
3 1/2	88mm	211 - 235mm	204 - 228mm	26-38mm	57
4	89mm	214 – 238mm	207 – 231mm	27-39mm	58
4 1/2	90mm	217 - 241mm	210 - 234mm	28-40mm	58
5	91mm	220 – 244mm	214 – 238mm	28-40mm	59
5 1/2	92mm	223 - 247mm	217 - 241mm	29-41mm	60
6	93mm	227 – 251mm	220 – 244mm	29-41mm	61
6 1/2	94mm	230 - 254mm	223 - 247mm	30-42mm	61
7	95mm	234 – 258mm	227 – 251mm	30-42mm	62
7 1/2	96mm	237 - 261mm	230 - 254mm	31-43mm	63
8	97mm	241 – 265mm	234 – 258mm	31-43mm	64
8 1/2	98mm	244 - 268mm	237 - 261mm	32-44mm	64
9	99mm	247 – 271mm	240 – 264mm	32-44mm	65
9 1/2	100mm	250 - 274mm	243 - 267mm	33-45mm	66
10	101mm	253 – 277mm	246 – 270mm	33-45mm	67

Very Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	88mm	211 – 235mm	202 – 226mm	26-38mm	54
1 1/2	89mm	214 - 238mm	205 - 229mm	27-39mm	55
2	90mm	217 – 241mm	208 – 232mm	27-39mm	56
2 1/2	91mm	220 - 244mm	211 - 235mm	28-40mm	56
3	92mm	223 – 247mm	215 – 239mm	28-40mm	57
3 1/2	93mm	226 - 250mm	218 - 242mm	29-41mm	58
4	94mm	229 – 253mm	221 – 245mm	29-41mm	59
4 1/2	95mm	232 - 256mm	224 - 248mm	30-42mm	59
5	96mm	235 – 259mm	228 – 252mm	30-42mm	60
5 1/2	97mm	238 - 262mm	231 - 255mm	31-43mm	61
6	98mm	242 – 266mm	234 – 258mm	31-43mm	62
6 1/2	99mm	245 - 269mm	237 -261mm	32-44mm	62
7	100mm	249 – 273mm	241 - 265mm	32-44mm	63
7 1/2	101mm	252 - 276mm	244 - 268mm	33-45mm	64
8	102mm	256 – 280mm	248 – 272mm	33-45mm	65
8 1/2	103mm	260 - 284mm	252 - 276mm	34-46mm	65
9	104mm	263 – 287mm	255 – 279mm	34-46mm	66
9 1/2	105mm	266 - 290mm	258 - 282mm	35-47mm	67
10	106mm	269 – 293mm	261 – 285mm	35-47mm	68

Super Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	93mm	227 – 251mm	214 – 238mm	27-39mm	55
1 1/2	94mm	230 - 254mm	217 - 241mm	28-40mm	56
2	95mm	233 – 257mm	220 – 244mm	28-40mm	57
2 1/2	96mm	236 - 260mm	223 - 247mm	29-41mm	57
3	97mm	239 – 263mm	226 – 250mm	29-41mm	58
3 1/2	98mm	242 - 266mm	229 - 253mm	30-42mm	59
4	99mm	245 – 269mm	233 – 257mm	30-42mm	60
4 1/2	100mm	248 - 272mm	236 - 260mm	31-43mm	60
5	101mm	251 – 275mm	240 – 264mm	31-43mm	61
5 1/2	102mm	254 - 278mm	243 - 267mm	32-44mm	62
6	103mm	258 – 282mm	246 – 270mm	32-44mm	63
6 1/2	104mm	261 - 285mm	249 - 273mm	33-45mm	63
7	105mm	264 – 288mm	252 – 276mm	33-45mm	64
7 1/2	106mm	267 - 291mm	255 - 279mm	34-46mm	65
8	107mm	271 – 295mm	258 – 282mm	34-46mm	66
8 1/2	108mm	274 - 298mm	261 - 285mm	35-47mm	66
9	109mm	277 – 302mm	264 – 288mm	35-47mm	67
9 1/2	110mm	281 - 305mm	268 - 292mm	36-48mm	68
10	111mm	284 – 309mm	271 – 295mm	36-48mm	69

The above charts refer to our standard lasts.

Pro-Fit Shallow Toe Last - reduce toe depth by 6mm

SIZING CHART | LADIES

Carewear Extra Deep Lasts (Extra 6mm Deeper)

Narrow & Extra Deep Wide & Extra Deep

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	78mm	189 - 213mm	182 - 205mm	27-39mm	52
1 1/2	79mm	192 - 216mm	185 - 209mm	27-39mm	53
2	80mm	195 - 219mm	188 - 212mm	28-40mm	54
2 1/2	81mm	198 - 222mm	191 - 215mm	28-40mm	54
3	82mm	201 - 225mm	194 - 218mm	29-41mm	55
3 1/2	83mm	205 - 229mm	198 - 222mm	29-41mm	56
4	84mm	208 - 232mm	210 - 225mm	30-42mm	57
4 1/2	85mm	211 - 235mm	204 - 228mm	31-43mm	57
5	86mm	215 - 239mm	208 - 232mm	31-43mm	58
5 1/2	87mm	218 - 242mm	211 - 235mm	32-44mm	59
6	88mm	222 - 246mm	214 - 238mm	32-44mm	60
6 1/2	89mm	225 - 249mm	217 - 241mm	33-45mm	60
7	90mm	228 - 252mm	220 - 244mm	33-45mm	61
7 1/2	91mm	231 - 255mm	223 - 247mm	34-46mm	62
8	92mm	234 - 258mm	226 - 250mm	34-46mm	63
8 1/2	93mm	237 - 261mm	229 - 253mm	35-47mm	63
9	94mm	241 - 265mm	233 - 257mm	35-47mm	64
9 1/2	95mm	244 - 268mm	236 - 260mm	36-48mm	65
10	96mm	247 - 271mm	239 - 263mm	36-48mm	66

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	83mm	206 - 230mm	200 - 224mm	30-42mm	53
1 1/2	84mm	210 - 234mm	203 - 227mm	30-42mm	54
2	85mm	213 - 237mm	206 - 230mm	31-43mm	55
2 1/2	86mm	216 - 240mm	210 - 234mm	31-43mm	55
3	87mm	220 - 244mm	213 - 237mm	32-44mm	56
3 1/2	88mm	223 - 247mm	216 - 240mm	32-44mm	57
4	89mm	226 - 250mm	219 - 243mm	33-46mm	58
4 1/2	90mm	229 - 253mm	222 - 246mm	34-45mm	58
5	91mm	232 - 256mm	226 - 250mm	34-45mm	59
5 1/2	92mm	235 - 259mm	229 - 253mm	35-47mm	60
6	93mm	239 - 263mm	232 - 256mm	35-47mm	61
6 1/2	94mm	242 - 266mm	235 - 259mm	36-48mm	61
7	95mm	246 - 270mm	239 - 263mm	36-48mm	62
7 1/2	96mm	249 - 273mm	242 - 266mm	37-49mm	63
8	97mm	253 - 277mm	246 - 270mm	37-49mm	64
8 1/2	98mm	256 - 280mm	249 - 273mm	38-50m	64
9	99mm	259 - 283mm	252 - 276mm	38-50mm	65
9 1/2	100mm	262 - 286mm	255 - 279mm	39-51mm	66
10	101mm	265 - 289mm	258 - 282mm	39-51mm	67

Very Wide & Extra Deep

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	88mm	223 - 247mm	214 - 238mm	32-44mm	54
1 1/2	89mm	226 - 250mm	217 - 241mm	33-45mm	55
2	90mm	229 - 253mm	220 - 244mm	33-45mm	56
2 1/2	91mm	232 - 256mm	223 - 247mm	34-46mm	56
3	92mm	235 - 259mm	227 - 251mm	34-46mm	57
3 1/2	93mm	238 - 262mm	230 - 254mm	35-47mm	58
4	94mm	241 - 265mm	233 - 257mm	35-47mm	59
4 1/2	95mm	244 - 268mm	236 - 260mm	36-48mm	59
5	96mm	247 - 271mm	240 - 264mm	36-48mm	60
5 1/2	97mm	250 - 274mm	243 - 267mm	37-49mm	61
6	98mm	254 - 278mm	246 - 270mm	37-49m	62
6 1/2	99mm	257 - 281mm	249 - 273mm	38-50mm	62
7	100mm	261 - 285mm	253 - 277mm	38-50mm	63
7 1/2	101mm	264 - 288mm	256 - 280mm	39-51mm	64
8	102mm	268 - 292mm	260 - 284mm	39-51mm	65
8 1/2	103mm	272 - 296mm	264 - 288mm	40-51mm	65
9	104mm	275 - 299mm	267 - 291mm	40-52mm	66
9 1/2	105mm	278 - 302mm	270 - 294mm	41-53mm	67
10	106mm	281 - 305mm	273 - 297mm	41-53mm	68

Super Wide & Extra Deep

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
1	93mm	239 - 263mm	226 - 250mm	33-45mm	55
1 1/2	94mm	242 - 266mm	229 - 253mm	34-46mm	56
2	95mm	245 - 269mm	232 - 256mm	34-46mm	57
2 1/2	96mm	248 - 272mm	235 - 259mm	35-47mm	57
3	97mm	251 - 275mm	238 - 262mm	35-47mm	58
3 1/2	98mm	254 - 278mm	241 - 265mm	36-48mm	59
4	99mm	257 - 281mm	245 - 269mm	36-48mm	60
4 1/2	100mm	260 - 284mm	248 - 272mm	37-49mm	60
5	101mm	263 - 287mm	252 - 276mm	37-49mm	61
5 1/2	102mm	266 - 290mm	255 - 279mm	38-50mm	62
6	103mm	270 - 294mm	258 - 282mm	38-50mm	63
6 1/2	104mm	273 - 297mm	261 - 285mm	39-51mm	63
7	105mm	276 - 300mm	264 - 288mm	39-51mm	64
7 1/2	106mm	279 - 303mm	267 - 291mm	40-52mm	65
8	107mm	283 - 307mm	270 - 294mm	40-52mm	66
8 1/2	108mm	286 - 310mm	273 - 297mm	41-53mm	66
9	109mm	280 - 313mm	276 - 300mm	41-53mm	67
9 1/2	110mm	293 - 317mm	280 - 304mm	42-54mm	68
10	111mm	296 - 320mm	283 - 307mm	42-54mm	69

SIZING CHART | MENS

Carewear Standard & Carewear Pro-Fit Lasts Narrow Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	85mm	205 – 229mm	200 – 223mm	23-35mm	58
4 1/2	86mm	208 - 232mm	203 - 227mm	23-35mm	59
5	87mm	211 – 235mm	206 – 230mm	24-36mm	60
5 1/2	88mm	214 - 238mm	209 - 233mm	24-36mm	60
6	89mm	217 – 241mm	212 – 236mm	25-37mm	61
6 1/2	90mm	220 - 224mm	214 - 238mm	25-37mm	62
7	91mm	224 – 248mm	217 – 241mm	25-37mm	63
7 1/2	92mm	227 - 251mm	220 - 244mm	25-37mm	64
8	93mm	231 – 255mm	223 – 247mm	26-38mm	64
8 1/2	94mm	234 - 258mm	226 - 250mm	27-39mm	65
9	95mm	238 – 262mm	229 – 253mm	27-39mm	66
9 1/2	96mm	241 - 265mm	232 - 256mm	27-39mm	67
10	97mm	245 – 269mm	236 – 260mm	28-40mm	68
10 1/2	98mm	248 - 272mm	239 - 263mm	28-40mm	68
11	99mm	251 – 275mm	242 – 266mm	29-41mm	69
11 1/2	100mm	254 - 278mm	245 - 269mm	29-41mm	70
12	101mm	258 – 282mm	248 – 272mm	29-41mm	71
12 1/2	102mm	261 - 285mm	250 - 274mm	30-42mm	72
13	103mm	265 – 289mm	253 – 277mm	30-42mm	72
13 1/2	104mm	268 - 292mm	255 - 279mm	31-43mm	73
14	105mm	272 – 296mm	258 – 282mm	31-43mm	74

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	91mm	219 – 243mm	214 – 238mm	26-38mm	59
4 1/2	92mm	222 - 246mm	217 - 241mm	26-38mm	60
5	93mm	225 – 249mm	220 – 244mm	27-39mm	61
5 1/2	94mm	228 - 252mm	223 - 247mm	27-39mm	61
6	95mm	231 – 255mm	226 – 250mm	28-40mm	62
6 1/2	96mm	234 - 258mm	228 - 252mm	28-40mm	63
7	97mm	238 – 262mm	231 – 255mm	28-40mm	64
7 1/2	98mm	241 - 265mm	234 - 258mm	29-41mm	65
8	99mm	245 – 269mm	237 – 261mm	29-41mm	65
8 1/2	100mm	248 - 272mm	239 - 262mm	30-42mm	66
9	101mm	251 – 275mm	242 – 266mm	30-42mm	67
9 1/2	102mm	254 - 278mm	245 - 269mm	30-42mm	68
10	103mm	258 – 282mm	248 – 272mm	31-43mm	69
10 1/2	104mm	261 - 285mm	250 - 274mm	31-43mm	69
11	105mm	265 – 289mm	253 – 277mm	32-44mm	70
11 1/2	106mm	268 - 292mm	255 - 279mm	32-44mm	71
12	107mm	272 – 296mm	258 – 282mm	32-44mm	72
12 1/2	108mm	275 - 299mm	261 - 285mm	33-45mm	73
13	109mm	278 – 302mm	264 – 288mm	33-45mm	73
13 1/2	110mm	280 - 304mm	266 - 290mm	34-46mm	74
14	111mm	283 – 307mm	269 – 293mm	34-46mm	75

Very Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	97mm	234 - 258mm	228 - 252mm	28-40mm	60
4 1/2	98mm	237 - 251mm	231 - 255mm	28-40mm	60
5	99mm	240 - 264mm	234 - 258mm	29-41mm	61
5 1/2	100mm	243 - 267mm	237 - 261mm	29-41mm	62
6	101mm	246 - 270mm	240 - 264mm	30-42mm	62
6 1/2	102mm	249 - 273mm	242 - 266mm	30-42mm	63
7	103mm	253 - 277mm	245 - 269mm	30-42mm	64
7 1/2	104mm	257 - 281mm	248 - 272mm	31-43mm	65
8	105mm	261 - 285mm	251 - 275mm	31-43mm	66
8 1/2	106mm	264 - 288mm	253 - 277mm	32-44mm	66
9	107mm	268 - 292mm	256 - 280mm	32-44mm	67
9 1/2	108mm	271 - 295mm	258 - 282mm	33-45mm	68
10	109mm	274 - 298mm	261 - 285mm	33-45mm	69
10 1/2	110mm	277 - 301mm	264 - 288mm	34-46mm	70
11	111mm	280 - 304mm	267 - 291mm	34-46mm	71
11 1/2	112mm	283 - 307mm	269 - 293mm	35-47mm	72
12	113mm	287 - 311mm	272 - 296mm	35-47mm	73
12 1/2	114mm	290 - 314mm	275 - 299mm	36-48mm	74
13	115mm	294 - 318mm	278 - 302mm	36-48mm	75
13 1/2	116mm	298 - 322mm	281 - 305mm	37-49mm	76
14	117mm	302 - 326mm	284 - 308mm	37-49mm	77

Super Wide

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	103mm	249 - 273mm	242 - 266mm	29-41mm	61
4 1/2	104mm	252 - 276mm	245 - 269mm	29-41mm	62
5	105mm	255 - 279mm	248 - 272mm	30-42mm	62
5 1/2	106mm	258 - 282mm	251 - 275mm	30-42mm	63
6	107mm	261 - 285mm	254 - 278mm	31-43mm	64
6 1/2	108mm	264 - 288mm	256 - 280mm	31-43mm	65
7	109mm	268 - 292mm	259 - 283mm	31-43mm	66
7 1/2	110mm	272 - 296mm	262 - 286mm	32-44mm	67
8	111mm	276 - 300mm	265 - 289mm	32-44mm	68
8 1/2	112mm	280 - 304mm	267 - 291mm	33-45mm	68
9	113mm	284 - 308mm	270 - 294mm	33-45mm	69
9 1/2	114mm	287 - 311mm	273 - 297mm	34-46mm	70
10	115mm	290 - 314mm	276 - 300mm	34-46mm	71
10 1/2	116mm	293 - 317mm	278 - 302mm	35-47mm	72
11	117mm	296 - 320mm	281 - 305mm	36-48mm	73
11 1/2	118mm	299 - 323mm	283 - 307mm	36-48mm	74
12	119mm	302 - 326mm	286 - 310mm	37-49mm	74
12 1/2	120mm	305 - 329mm	288 - 312mm	38-50mm	75
13	121mm	308 - 332mm	290 - 314mm	38-50mm	76
13 1/2	122mm	311 - 335mm	293 - 317mm	39-51mm	77
14	123mm	314 - 338mm	295 - 319mm	39-51mm	78

The above charts refer to our standard lasts. Pro-Fit Shallow Toe Last - reduce toe depth by 6mm

SIZING CHART | MENS

Carewear Extra Deep Lasts (Extra 6mm Deeper)

Narrow & Extra Deep Wide & Extra Deep

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	85mm	217 - 241mm	212 - 236mm	29-41mm	58
4 1/2	86mm	220 - 244mm	215 - 239mm	29-41mm	59
5	87mm	223 - 247mm	218 - 242mm	30-42mm	60
5 1/2	88mm	226 - 250mm	221 - 245mm	30-42mm	60
6	89mm	229 - 253mm	224 - 248mm	31-43mm	61
6 1/2	90mm	232 - 256mm	226 - 250mm	31-43mm	62
7	91mm	236 - 260mm	229 - 253mm	31-43mm	63
7 1/2	92mm	239 - 263mm	232 - 255mm	32-44mm	64
8	93mm	243 - 267mm	235 - 259mm 32-44mm		64
8 1/2	94mm	246 - 270mm	238 - 262mm	33-45mm	65
9	95mm	250 - 274mm	241 - 265mm	33-45mm	66
9 1/2	96mm	253 - 277mm	244 - 268mm	33-45mm	67
10	97mm	257 - 277mm	248 - 272mm	34-46mm	68
10 1/2	98mm	260 - 284mm	251 - 275mm	34-46mm	68
11	99mm	263 - 287mm	254 - 278mm	35-47mm	69
11 1/2	100mm	266 - 290mm	257 - 281mm	35-47mm	70
12	101mm	270 - 294mm	260 - 284mm 35-47m		71
12 1/2	102mm	273 - 297mm	262 - 285mm	36-48mm	72
13	103mm	277 - 301mm	265 - 289mm	36-48mm	72
13 1/2	104mm	280 - 304mm	267 - 291mm 37-49mn		73
14	105mm	284 - 308mm	270 - 294mm	37-49mm	74

Size	Joint Width	Joint Circumference	Instep Circumference	Toe Height Min/Max	Heel Width
4	91mm	231 - 255mm	226 - 250mm	32-44mm	59
4 1/2	92mm	234 - 258mm	229 - 253mm	32-44mm	60
5	93mm	237 - 261mm	232 - 256mm	33-45mm	61
5 1/2	94mm	240 - 264mm	235 - 259mm	33-45mm	61
6	95mm	243 - 267mm	238 - 262mm	34-46mm	62
6 1/2	96mm	246 - 270mm	240 - 264mm	34-46mm	63
7	97mm	250 - 274mm	243 - 267mm	34-46mm	64
7 1/2	98mm	253 - 277mm	246 - 270mm	35-47mm	65
8	99mm	257 - 281mm	249 - 273mm	35-47mm	65
8 1/2	100mm	260 - 284mm	251 - 275mm	36-48mm	66
9	101mm	263 - 287mm	254 - 278mm	36-48mm	67
9 1/2	102mm	266 - 290mm	257 - 281mm	36-48mm	68
10	103mm	270 - 294mm	260 - 284mm	37-49mm	69
10 1/2	104mm	273 - 297mm	262 - 285mm	37-49mm	69
11	105mm	277 - 301mm	265 - 289mm	38-50mm	70
11 1/2	106mm	280 - 304mm	267 - 291mm	38-50mm	71
12	107mm	284 - 308mm	270 - 294mm	38-50mm	72
12 1/2	108mm	287 - 311mm	273 - 297mm	39-51mm	73
13	109mm	290 - 314mm	276 - 300mm	39-51mm	73
13 1/2	110mm	293 - 317mm	278 - 302mm	40-52mm	74
14	111mm	296 - 320mm	281 - 305mm	40-52mm	75

Very Wide & Extra Deep

Super Wide & Extra Deep


Size	Joint Width			Toe Height Min/Max	Heel Width
4	97mm	246 - 270m	240 - 264mm	34-46mm	60
4 1/2	98mm	249 - 273mm	243 - 267mm	34-46mm	60
5	99mm	252 - 276mm	246 - 270mm	35-47mm	61
5 1/2	100mm	255 - 279mm	249 - 273mm	35-47mm	62
6	101mm	258 - 282mm	252 - 276mm	36-48mm	62
6 1/2	102mm	261 285mm	254 - 278mm	36-48mm	63
7	103mm	265 - 289mm	257 - 281mm	36-48mm	64
7 1/2	104mm	269 - 293mm	260 - 284mm	37-49mm	65
8	105mm	273 - 297mm	263 - 287mm	37-49mm	66
8 1/2	106mm	276 - 300mm	265 - 289mm	38-50mm	66
9	107mm	280 - 304mm	268 - 292mm	38-50mm	67
9 1/2	108mm	283 - 307mm	270 - 294mm	39-51mm	68
10	109mm	286 - 310mm	273 - 297mm	39-51mm	69
10 1/2	110mm	289 - 313mm	276 - 300mm	40-52mm	70
11	111mm	292 - 316mm	279 - 303mm	40-52mm	71
11 1/2	112mm	295 - 319mm	281 - 305mm	41-53mm	72
12	113mm	299 - 323mm	284 - 308mm	41-53mm	73
12 1/2	114mm	303 - 327mm	287 - 311mm	42-54mm	74
13	115mm	307 - 331mm	290 - 314mm	42-54mm	75
13 1/2	116mm	311 - 335mm	293 - 317mm	43-55mm	76
14	117mm	314 - 338mm	296 - 320mm	43-55m	77

Size	Joint Width	Joint Instep Circumference Circumference		Toe Height Min/Max	Heel Width
4	103mm	261 - 285mm	254 - 278mm	35-47mm	61
4 1/2	104mm	264 - 288mm	257 - 281mm	35-47mm	62
5	105mm	270 - 294mm	260 - 284mm	36-48mm	62
5 1/2	106mm	270 - 294mm	263 - 287mm	36-48mm	63
6	107mm	273 - 297mm	266 - 290mm	37-49mm	64
6 1/2	108mm	276 - 300mm	268 - 292mm	37-49mm	65
7	109mm	280 - 304mm	271 - 295mm	37-49mm	66
7 1/2	110mm	284 - 308mm	274 - 298mm	38-50mm	67
8	111mm	288 - 312mm	277 - 301mm	38-50mm	68
8 1/2	112mm	292 - 316mm	279 - 303mm	39-51mm	68
9	113mm	296 - 320mm	282 - 306mm	39-51mm	69
9 1/2	114mm	299 - 223mm	285 - 309mm	40-52mm	70
10	115mm	302 - 326mm	288 - 312mm	40-52mm	71
10 1/2	116mm	305 - 329mm	290 - 314mm	41-53mm	72
11	117mm	308 332mm	293 - 317mm	42-54mm	73
11 1/2	118mm	311 - 335mm	295 - 319mm	42-54mm	74
12	119mm	314 - 338mm	298 - 322mm	43-55mm	74
12 1/2	120mm	317 - 341mm	300 - 324mm	44-56mm	75
13	121mm	320 - 344mm	303 - 327mm	44-56mm	76
13 1/2	122mm	323 - 347mm	305 - 329mm	45-57mm	77
14	123mm	326 - 350mm	308 - 332mm	45-57mm	78

Our Carewear Last

Our most popular last, the Carewear last provides a comprehensive range of sizes, width fittings and depths.


Available as: Standard, Pro-Fit and Extra Deep



Carewear Standard last

Available in four widths. Narrow, Wide, Very Wide and Super Wide. Women's size 1 to 10. Men's size 4 to 14. Including half sizes.

Allowance for up to 12mm of cushioning.




Carewear Pro-Fit last

The Pro-Fit last features a reduced toe depth (6mm) for an improved cosmetic appearance and a tighter back curve and top line to enhance the fit at the back of the shoe.

The Pro Fit Last retains the Standard lasts joint and instep dimensions, giving patients the width and depth where they need it the most. Available in four widths. Narrow, Wide, Very Wide and Super Wide. Women's size 1 to 10. Men's size 4 to 14, including half sizes.

Allowance for up to 12mm of cushioning.



Carewear Extra Deep last

Retaining the same toe shape and profile of the Standard last, the Extra Deep version features an additional 6mm depth. Available in four widths, Narrow, Wide, Very Wide and Super Wide. Women's Sizes 1-10 and men's Sizes 4 to 14, including half sizes.

Suitable for more extreme foot problems such as oedema and extra high retracted toes, or where extra cushioning is required with the use of moulded total contact insoles.

Allowance for up to 18mm of cushioning.

KEN HALL | OPTIONS

Upper Materials











Lining Materials Leather Plush









COLOUR CHART

Leather | Ultra soft leather when extreme comfort is required



Nubuck | Soft with a milled finish and tight nap



Suede | Classic nap finish





Green



Waxy Robust but with a soft handle for orthopaedic use. Additional colours available on request.









MODULAR ORDER FORM

DATE	HOSPITAL	
TRUST	CLINICIAN	
ORDER NUMBER	CLINICIAN EMAIL	
PATIENT NAME	CLINICIAN PHONE	

LAST DETAILS							
SIZE & MEASUREMENTS	LEFT	RIGHT	OTHER LAST DETAILS				
SIZE:							
JOINT CIRCUMFERENCE:							
INSTEP CIRCUMFERENCE:							
JOINT WIDTH:							
HEEL WIDTH:			INSOLE / TCI SPECIFICATIONS				
LAST (PLEASE TICK)	LEFT	RIGHT	LEFT				
CAREWEAR	TICK	TICK					
CAREWEAR EXTRA DEEP	TICK	TICK					
CAREWEAR PROFIT	TICK	TICK					
WIDTH (PLEASE TICK)	LEFT	RIGHT	RIGHT				
NARROW	TICK	TICK					
WIDE	TICK	TICK					
VERY WIDE	TICK	TICK					
SUPER WIDE	TICK	TICK					

UPPER DETAILS								
STYLE			OTHER DETAILS					
COLOUR								
FASTENING (PLEASE TICK)	LACE	VELCRO						
(FLEASE HCK)	TICK	TICK						
UPPER MATERIAL	SOFTEE	CALF		NUBUCK	SUEDE	AS CAT	ALOGUE	
(PLEASE TICK)	TICK	TICK		TICK TICK		TICK		
LINING MATERIAL (PLEASE TICK)	PLUSH	STRAW		PIG	SHEARLING			
(FLEASE HOR)	TICK	TICK		TICK TICK		K		
TOE STIFFENER	WING	FULL		MEMORY	STRONG TOE PUFF STRONG		STRONG WING	
(PLEASE TICK)	TICK	TICK		TICK	TICK		TICK	
LEFT HEEL STIFFENER			RIGHT	HEEL STIFFENER				

SOLE DETAILS								
SOLE DESIGN (PLEASE TICK)	SEPARATE	WEDGE	SCALLOP	EVA COMMANDO		TRAINER		
(FLEASE HOK)	TICK	TICK	TICK	TICK		TICK		
SOLE SPECIFICATIONS	HARD COMMANDO		TICK	BLAKED	TICK	WELTED	TICK	
	•			•	-			



01536 522468 info@kenhall.co.uk www.kenhall.co.uk

Ken Hall Ltd Newman Street Kettering NN16 0TG

- in /ken-hall
- /KenHallLtd
- @kenhallfootwear



